



LAMPE OPERATOIRE

Série OL-200

Lampes opératoires – scialytique HALOGENE 100.000 LUX mono ampoule



Style No.	Intensity	Color Temp.	Illuminate Field	Vertical Adjustment	Lamp-head Diameter	Bulb
	Lux	k	φ mm	mm	φ mm	v/w
OL-201	100,000	4,000± 500	100-200	900	700	24/150 halogen
OL-202	Main 100,000	4,000± 500	100-200	900	700	24/150 halogen
	Aux.100,000	4,000± 500	100-200	900	700	24/150 halogen
OL-203	100,000	4,000± 500	100-200	900	600	24/150 halogen
OL-204	Main 100,000	4,000± 500	100-200	900	600	24/150 halogen
	Aux.100,000	4,000± 500	100-200	900	600	24/150 halogen

