



MONITEUR MULTIPARAMETRES POLYVALENT

Patient Monitor

Appareil de nouvelle génération, polyvalent et complet
Existe en de nombreuses versions (Réanimation, Anesthésie,...)
Prévu pour adulte et nouveau-né

Standard Feature:

12.1 inch color TFT display,
4-channel ECG, ST, Arrhythmia, RESP, SPO2, NIBP, HR, PR,
TEMP
Drug dose calculation,
Pacemaker detection,
Apnea alarm,
CNS interface,
Rechargeable internal battery

Standard accessories:

1x ECG 5-lead wire & cable, 1x pack Electrodes, 1x SPO2
probe, 1x NIBP cuff, T° probe

Optional:

2xIBP, CO₂, CO, Multi-gas, Low perfusion (SAC), 3-channel built-in thermal recorder

Specification:

12.1 inches Color TFT Display;
Big font display;
Drug Dose Calculation,
Maximum 8-channel Waveform Display (4x ECG
Channels)
600-set Blood Pressure Data Storage and Recall;
120 hours Trend information;
ST Segment, Arrhythmia Analysis
Pacemaker Detection and Apnea Alarm



12 Minutes Real Time Waveform History Review
Hemodynamic calculation
Manual/Auto/Stat Measurement of NIBP
Applications from Neonate to Adult;
External display function
Networking Capability
AC/DC Operation
Build-in battery