



NEBULISEUR AEROSOL MIKO

MIKO is an electrical feed compressor with aerosol therapy nebulizer. The instrument is designed for easy transport and handling and it's recommended for atomising antibiotics and bronchodilator drugs. The high thermal insulated plastic body complies with European Safety Standards. The instrument is designed for not continuous use 20 min. ON / 40 min. OFF.



TECHNICAL SPECIFICATIONS	
Power Feeding (other voltages and plugs available on request)	230V - 50Hz European 2-pin plug / spina Italiana a 2 poli
Power Consumption	170 VA
Fuse	F1 x 1.6 A 250 V
Max Pressure	2.2 bar - 32 psi - 220 kPa
Operating Pressure	0.95 bar - 13.8 psi - 95 kPa
Max Air Flow	13 l/min
Operating Air Flow	4,8 l/min
Duty Cycle	20 min ON - 40 min OFF
Noise Level (EN 13544-1 normative)	Approx. 56 dBA
Weight	3.08 Lbs - 1.45 kg
Size	in 5,83" x 4,49" x 6,73" - mm 148 x 114 x 171
Years of Warranty	3
Shipping carton	6
Manufacturing	100% in Italy
AEROSOL SPECIFICATIONS	
Neb-Rate with 4 ml of solution	0,35 ml/min x 0,9% NaCl
Mass Median Aerodynamic Diameter (MMAD) measured in compliance with EN 13544-1	2,44 µm
Respirable Fraction (<5µm)	> 80%