



VERTICAL AUTOCLAVE INOX

Vertical type stainless steel pressure steam sterilizer with Hand wheel type of quick-open door

Model: of 40, 50, 75, 100, 150 and 200 litres
Quality: Stainless steel sterilizing chamber

Product Description:

- Hand wheel type of quick-open door structure
- Computer control auto recycle sterilization
- LDC screen indicate working state
- Pressure safety lock
- Temperature scope 50-126°C
- Time scope 0-99h
- Overpressure auto-discharging 0.145-0.165MPa
- Automatic shut down with beep reminding after sterilization
- Stainless steel sterilizing chamber



| Type | Volume | Power source | Consumed power | Dimension of sterilizing chamber | Weight | Dimension |
|------|--------|--------------|----------------|----------------------------------|--------|--------------|
| 40 | 40L | 220V | 3.5Kw | (Φ)350x450 | 80kg | 600x580x1250 |
| 50 | 50L | 220V | 3.5Kw | (Φ)350x520 | 80kg | 600x580x1250 |
| 75 | 75L | 220 V/380V | 3.5Kw/4.5Kw | (Φ)388x650 | 90 kg | 600x600x1400 |
| 100 | 100L | 380 V | 4.5 Kw | (Φ)500×550 | 100 kg | 630×630×1300 |
| 150 | 150L | 380 V | 6.0 Kw | (Φ)500×760 | 110 kg | 630×630×1450 |
| 200 | 200L | 380 V | 6.0 Kw | (Φ)600×700 | 120 kg | 730×730×1390 |