



# VERTICAL AUTOCLAVE INOX

## Vertical type stainless steel pressure steam sterilizer with Hand wheel type of quick-open door

Model: of 40, 50, 75, 100, 150 and 200 litres  
Quality: Stainless steel sterilizing chamber

### Product Description:

- Hand wheel type of quick-open door structure
- Computer control auto recycle sterilization
- LDC screen indicate working state
- Pressure safety lock
- Temperature scope 50-126°C
- Time scope 0-99h
- Overpressure auto-discharging 0.145-0.165MPa
- Automatic shut down with beep reminding after sterilization
- Stainless steel sterilizing chamber



Type	Volume	Power source	Consumed power	Dimension of sterilizing chamber	Weight	Dimension
40	40L	220V	3.5Kw	(Φ)350x450	80kg	600x580x1250
50	50L	220V	3.5Kw	(Φ)350x520	80kg	600x580x1250
75	75L	220 V/380V	3.5Kw/4.5Kw	(Φ)388x650	90 kg	600x600x1400
100	100L	380 V	4.5 Kw	(Φ)500x550	100 kg	630x630x1300
150	150L	380 V	6.0 Kw	(Φ)500x760	110 kg	630x630x1450
200	200L	380 V	6.0 Kw	(Φ)600x700	120 kg	730x730x1390