



BRONCHOSCOPE XZ-10WP

Broncho - fiberscope XZ-10 WP CE

The Bronchofiberscope is intended for observation, diagnosis and treatment of diseases of trachesbronchial tree, including biopsy, brush cytology, demonstration and photography. With the attachment of the suction adapter, this instrument can also be used for suction, injection and H.P.D., etc...

XZ-10 is the newly launched model of broncho-fiberscope with super thinner insertion tube:3.6mm and water-proofed.

Thanks for it's special technical spec., this model is suitable for the examination of pathological changes in sub-peripheral areas of the lungs and it also be used to make static and dynamic records, teaching by demonstration system (OES)



First Class Optics:

The 85° angle of view improving light system enhances orientation and helps eliminate blind spots. The images are large and clear and the same quality image can be gained when using TV system.

New Shape Control Section:

The shape of the control section has been improved to fit the hand more comfortably while permitting continuous uninterrupted movement of the fingers. The forceps inlet removes to the lower part of the control section, with helps the physician greatly to carry out the biopsy.

Special Insertion Tube Flexibility:

Ease to introduction is enhanced by the combination of an atraumatic tube with specially designed flexibility characteristics.

Good Compatibility:

The fiberscope has a compatibility with Olympus cold light source

Specifications:

Angle of view:	85°
Depth of view field:	3-50 mm
Outer diameter of insertion tube:	3.6 mm
Outer diameter of distal end	3.6 mm
Bending capability Up / Down:	130° / 90°
Inner diameter of biopsy channel	2.0 mm
Working length:	550 mm
Total length:	790 mm



Supplied with:

1	Sampling forceps for biopsy	1	Mouthpiece
1	Cytological brush	1	Forceps seal valve
1	Cleaning brush	1	Injection tube
1	Suction tube		