



# CENTRIFUGEUSE EASY

## ■ EASY CENTRIFUGE C-12B

- Suitable for laboratory and clinical tests.
- Adjustable speed design.
- Can choose with 6ml or 10 ml adapter
- Power pilot light & fuse safety device.
- Manual cover lock.



Model	Speed (rpm)	Max. RCF (x g)	Max. Capacity (ml)	Weight (kg)	Size (cm) W x D x H
C-12B	0 ~ 3.000	900	4x10 or 4x6	3	25 x 25 x 25

AC 110 / 220V, 50 / 60Hz