



## CHOLEDOCHOFIBERSCOPE

### Choledo - fiberscope CHS-1

Model CHS-1 is one of the principal instruments to make diagnosis and treatment of cystic duct disease. It is also an instrument of originality for treatment of post cystic duct and hepatic residual store. It helps patient to avoid suffering from surgical operation once more.

Ease of introduction is enhanced by the combination of an atraumatic insertion tube with specially designed flexibility characteristics.

There are fluid filling, ware ring and suction valves built in at the operating handle . With the advantage of large aperture diameter of forceps channel where even if instrument is handling, fluid filling, watering and suction may be proceeded as well.

#### First Class Optics:

The 85° angle of view improving light system enhances orientation and helps eliminate blind spots. The images are large and clear and the same quality image can be gained when using TV system.

#### Good Compatibility:

The fiberscope has a compatibility with Olympus cold light source

#### Specifications:

<b>Angle of view:</b>	90°
<b>Depth of view field:</b>	3-50 mm
<b>Outer diameter of insertion tube:</b>	6 mm
<b>Outer diameter of distal end</b>	5.8 mm
<b>Bending capability Up/Down:</b>	130°
<b>Bending capability Right/Left:</b>	100°
<b>Inner diameter of biopsy channel</b>	2.5 mm
<b>Working length:</b>	330 mm
<b>Total length:</b>	650 mm

