



COLONO FIBROSCOPE

Colonoscope CS-10 (CE)

First Class Optics:

The 100° angle of view improving light system enhances orientation and helps eliminate blind spots. The images are large and clear and the same quality image can be gained when using TV system.

Special Insertion Tube Flexibility:

The structure of the angle control is designed to fit the hand more ingenious. With the advanced technology, the insertion tube has specially flexibility characteristics. It can be inserted easily, decreasing the tiredness of the doctors and the pain of the patient.

Good Compatibility:

The fiberscope has a compatibility with Olympus cold light source (OES)



Specifications:

Field of view:	100°
Observation range:	5-100 mm
Viewing direction:	Forward
Distal end diameter:	12.9 mm
Forceps channel diameter:	3.2 mm
Bending capability Up/Down:	180°/180°
Bending capability Right/Left:	160°/160°
Working length:	1,3500 mm
Total length:	1,670 mm

Supplied with:

- 2 Sampling forceps for biopsy
- 1 Cytological brush
- 1 Cleaning brush
- 1 Infusion tube
- 1 Water container
- 5 Protective cap for forceps channel
- 8 Shaped seal O Ring

