



# COLONO FIBROSCOPE

## Colonoscope CS-120 (CE)

### First Class Optics:

The 100° angle of view improving light system enhances orientation and helps eliminate blind spots. The images are large and clear and the same quality image can be gained when using TV system.

### Special Insertion Tube Flexibility:

The structure of the angle control is designed to fit the hand more ingenious. With the advanced technology, the insertion tube has specially flexibility characteristics. It can be inserted easily, decreasing the tiredness of the doctors and the pain of the patient.



### Good Compatibility:

The fiberscope has a compatibility with Olympus cold light source (OES)

### Specifications:

<b>Field of view:</b>	100°
<b>Observation range:</b>	5-100 mm
<b>Viewing direction:</b>	Forward
<b>Distal end diameter:</b>	13.5 mm
<b>Forceps channel diameter:</b>	3.2 mm
<b>Bending capability Up/Down:</b>	180°/180°
<b>Bending capability Right/Left:</b>	160°/160°
<b>Working length:</b>	1,200 mm
<b>Total length:</b>	1,500 mm

### Supplied with:

- 2 Sampling forceps for biopsy
- 1 Cytological brush
- 1 Cleaning brush
- 1 Infusion tube
- 1 Water container
- 5 Protective cap for forceps channel
- 8 Shaped seal O Ring

