



COLPOSCOPE COLD LIGHT

No Compromise On Quality

Successful diagnosis begins with the colposcope. This colposcope, with its modern, contemporary design, has been specially developed for daily diagnosis purposes.

Despite its attractive price, it meets our high standards of quality, while not compromising accessories.

A comprehensive range of accessories is available for this handy, flexible instrument. Whether mounted to a stand, the high-power, fibre optic, cold light source combined with single magnification or 3-step changer, straight or inclined binocular tubes, interchangeable lenses and eyepieces, the T-handle which can be mounted at either side, provide an exceptionally high degree of flexibility in day-to-day use.

An optional video camera can be integrated enabling co-observation as well as documentation via video tape recording or computerization.



Technical detail :

Rolling base	600 mm height
Binocular straight tube	f = 125 mm
Binocular inclined tube 45°	f = 125 mm
Objective	f = 300 mm
Eyepieces (WF)	WF 20 x V (on choice WF 16 x V)
Green filter	swing in green filter
Cold light illumination	12V/75W continuously adjustable (110V and 220V)
T-handle	can be mounted left of right side

