



## BALANCE FOR LABORATORIES

### Compact balance

The balance has an AUTO-OFF function after 3 min to preserve the battery.

The Socket for the mains adapter as standard.

Large LCD display with digit height 16 mm.

Weight in function simplifies mixing.

Total balance W X D X H 170 X 240 X 38 mm.

Net weight of the balance approx. 600g.



Different available models:

Model	Max weight	Accuracy	Φ (mm)
A	200 g	0.1 g	150
D	400 g	0.1 g	125
B	2 kg	1 g	150
C	5 kg	5 g	150