



## CYSTOSCOPE

### Cystoscope *CYS-1*

Model *CYS-1* is one of the principal instruments to make diagnosis and treatment of pathologic changes in bladder. It helps patient to avoid suffering from surgical operation once more.

Ease of introduction is enhanced by the combination of an atraumatic insertion tube with specially designed flexibility characteristics.

There are fluid filling, ware ring and suction valves built in at the operating handle . With the advantage of large aperture diameter of forceps channel where even if instrument is handling, fluid filling, watering and suction may be proceeded as well.

#### **First Class Optics:**

The 75° angle of view improving light system enhances orientation and helps eliminate blind spots. The images are large and clear and the same quality image can be gained when using TV system.

#### **Good Compatibility:**

The fiberscope has a compatibility with Olympus cold light source

#### **Specifications:**

<b>Angle of view:</b>	75°
<b>Outer diameter of insertion tube:</b>	6,3 mm
<b>Bending capability Up/Down:</b>	180°
<b>Bending capability Right/Left:</b>	90°
<b>Inner diameter of biopsy channel</b>	2 mm
<b>Working length:</b>	380 mm

