



FCETAL DOPPLER PD1+

Description :



Spécifications :

Ultrasound

Frequency	2 MHz continuous wave
Transducer	2 crystal narrow beam
Output power	<50mW/mm ² SATA
Audio response	300Hz - 1KHz
Fetal heart rate	Multipoint real time autocorralator

Unit Controls

Keys	1 Key (for unit on/off)
Controls	Rotary volume
Indicators	3 digital FHR LCD display FHR pulse icon Battery Low icon

Power Supply

Battery	MNI604 (PP3) 9v Alkaline Manganese
Expected battery life	>6 hours

Outputs

Headset	Audio output to optional headset
Serial	RS232 interface to optional UltraTrace 2 PC software

Enclosure

Material	ABS / Noryl
Size	150mm by 75MM
Weight	350g typically including transducer

Safety

Classification	Class 1 Type B - IEC 60601-1
----------------	------------------------------