



FETAL DOPPLER PD1

Description :



Spécifications :

Ultrasound

Frequency	2 MHz continuous wave
Transducer	2 crystal narrow beam
Output power	<50mW/mm ² SATA
Audio response	300Hz - 1KHz

Unit Controls

Keys	1 Key (for unit on/off)
Controls	Rotary volume
Indicators	Yellow battery low LED

Power Supply

Battery	MNI604 (PP3) 9v Alkaline Manganese
Expected battery life	>9 hours

Outputs

Headset	Audio output to optional headset
---------	----------------------------------

Enclosure

Material	ABS / Noryl
Size	150mm by 75mm
Weight	350g typically including transducer

Safety

Classification	Class 1 Type B - IEC 60601-1
----------------	------------------------------