



LARYNGO FIBROSCOPE

Laryngoscope XH-3 (CE) WATERPROOF

First Class Optics:

The 85° angle of view, 24 x eyepiece and improving light system enhances orientation and helps eliminate blind spots. The images are large and clear and the same quality image can be gained when using TV system.



Special Insertion Tube Flexibility:

The structure of the angle control is designed to fit the hand more ingenious. With the advanced technology, the insertion tube has specially flexibility characteristics. It can be inserted easily, decreasing the tiredness of the doctors and the pain of the patient.

Good Compatibility:

The fiberscope has a compatibility with Olympus cold light source (OES)

Specifications:

Field of view:	85°
Observation range:	3-50 mm
Viewing direction:	Forward
Distal end diameter:	5.8 mm
Forceps channel diameter:	2.0 mm
Bending capability Up/Down:	130°/90°
Working length:	410 mm
Total length:	650 mm
Purge of tube whole:	yes

Supplied with:

- 1 Forceps for biopsy
- 1 Cytological brush
- 1 Cleaning brush
- 1 Suction tube
- 1 Mouthpiece
- 1 Forceps seal valve
- 1 Injection tube
- 1 Sterilized connector
- 1 Leak measure

