



COLPOSCOPE Leisegan 1D / 1DS

Standard Colposcopes

Excellent optics for highest-quality images
Fatigueless, simple and time-saving operation
Bright, continuously variable illumination
Large selection of accessories

Technical details:

Visual magnifications:

7.5 / 15 / 30 x (Model 1D)

3.75 / 7.5 / 15 x (Model 1DS)

Viewing field, diameter at:

3.75 x magnification: 90 mm

7.5 x magnification: 45 mm

15 x magnification: 22.5 mm

30 x magnification: 11 mm


Working distance : 300 mm

Settings : Focus : 40 mm
Height control: 80 mm

Illumination: 12V / 20W

Filter: green

Compensation by dioptré: +7 to -7 diopters



Modell 1D / Model 1D

Available stands



Balance-o-matic mount

/ Tilt stand

/ Upright stand