




PCA SYRINGE PUMP A5805


Ideal for :

- Pain Control Therapy
- Chemotherapy in Oncology
- Patient transfer or rescue
- and many more...



Specifications:

Power Supply	DC 4.5V (3 AA batteries)
Applicable syringes	2ml - 12ml disposable syringes
Flow Rate	1 mm/hr - 99 mm/hr (equals to 0.2 ml/hr - 20 ml/hr)
Syringe Accuracy	less than $\pm 2\%$
Alarms	Audible and Visual Alarms of "COMPLETE / OCCLUSION" and "LOW BATT"
Included Accessories	<ul style="list-style-type: none"> a) 4 AA batteries b) PCA bolus cable c) Syringe protecting cover c) Leather carrying belt
P.C.A. Function	Each PCA (Patient-Controlled Analgesia) dosage is 1/4 of preset delivery dosage. The interval between two P.C.A injections is no less than 30 minutes
Dimensions	L x B x H: 166x60x30 mm
Weight	250 g / 0.55 lb
Applications	<p>Can be used in variety of medicine therapies, great for adjustable long term infusion of medication. Indicative application fields are :</p> <ul style="list-style-type: none"> - Pain Control Therapy , - Patient transfer or rescue, Immunoglobulin Infusion, Heparin Administration , Chemotherapy in Oncology , Coronary & Neonatal Care , Treatment in Thalassaemia , Gynecology & Obstetrics



Portable/light in weight



Convenient carrying leather case & PCA bolus cable