



PCA AMBULATORY INFUSION PUMP A5806

This ambulatory infusion pump delivers a constant, controlled infusion rate of local anesthetic directly into the operative site for a relatively pain-free recovery. Patient can push PCA (Patient Controlled Analgesia) button to give themselves the extra dosage to address the breakthrough pain.

Ideal Solution for

- Labor Epidurals
- Regional Anesthesia
- Wound Management
- Chronic and Acute Pain
- Post-Operative Pain Management



Simple, Convenient and Affordable

- Simple Operation
- Simply Dial to Set Infusion Rate in Seconds
- Complete Portable (carrying belt included)
- Small and Light in Weight
- Uses AA Batteries (included)
- For Both Hospital and Home-care Use
- Pump is affordable for Disposable Use

Ambient temperature	5°C - 40°C
Power supply	AA battery x 2
Medication Reservoirs	100 ml or 260 ml
Flow rate option	Purge; 1, 2, 3, 4, 6, 8, 10 ml/hr
Alarms	Audible and Visual Alarms of "COMPLETE / ERROR" and "LOW BATT"
P.C.A. Function	Each PCA (Patient-Controlled Analgesia) dosage is 1/4 of preset delivery dosage. The interval between two P.C.A injections is no less than 30 minutes
Included Accessories w/ each Reservoir	a) Disposable Reservoir Box (sterilized package) b) AA Battery x 2 c) Disposable Extension Tubing
Dimensions	L x B x H : 140x95x35mm
Weight	200g / 0.44 lb (one pump with one 100ml reservoir)