



PHARMACEUTICAL REFRIGERATOR

Features:

1. Microprocessor temperature control system with digital the temperature display precisely controls the temperature inside the cabinet, the temperature the cabinet is clear at a glance.
2. Quality steel wire shelves soaked with lacquer provides convenient access to the stored goods and easy cleaning.
3. Quality carbon steel casing coated with imported paint by advanced coating technology provides surface color and lustre match well, collision resistance and corrosion resistance.
4. 90mm extra environment friendly insulation reduces effectively loss of cold air and provides energy saving and prolongs service life of the refrigerator.
5. Interior cabinet is made of coated with plastic, galvanized steel sheet.
6. Automatic defrost function, more convenient.
7. Use of refrigerant of R134a, a true green product.
8. Perfect audible/visual alarm system.



Specifications:

Model	Y-100	Y-200	Y-300	Y-660
Type	Upright			
Capacity	100 L	200 L	300 L	660 L
Door	1	1	1	2
Climate class	N	N	N	N
Protective class	I/B	I/B	I/B	I/B
Temperature	0°C to 10°C adjustable			
Temperature control precision	±1	±1	±1	±1
Mode of alarm	Audible/Visual alarm			
Power supply	220V/50Hz	220V/50Hz	220V/50Hz	220V/50Hz
Defrost mode	Automatic			
Rated input power	150W	300W	300W	550W
Power consumption	1.80 kwh/24h	2.80 kwh/24h	3.20 kwh/24h	4.00 kwh/24h
Refrigerant	R134a	R134a	R134a	R134a
Net weight	55 Kg	95 Kg	130 Kg	175 Kg
Gross weight	63 Kg	110 Kg	152 Kg	215 Kg
External dimension (mm)	540×615×1085	540×615×1725	610×680×1885	1100×650×2050
Packing dimension (mm)	600×680×1100	600×680×1740	666×750×1900	1150×720×2150
Loading quantity	114 pcs/40GP	57 pcs/40GP	51 pcs/40GP	30 pcs/40GP
Certificate	CE	CE	CE	CE