



MOBILE RX PL 102

PL102 Mobile X-ray



Power Output :	3KW
Commute Method	Full Wave Commutation
X-Ray Tube :	Fixed anode
Tube Voltage :	40~90kV
Tube Current :	50~90KV 15mA
	50~90KV 30mA Tube Current
	40~80KV 50mA
	50~90KV 15mA 0.1~6.3s
Exposure time	50~90KV 30mA 0.1~6.3s
	40~80KV 50mA 0.1~1.5s
mAs :	1.0~125 mAs (45 steps)
Power Supply :	220V±10% / 50Hz inner-resistance ≤1.0Ω
Operation Method :	Wire/Wireless control
Packing Details:	1240mm×750mm×1460mm
	N.W.150kgs G.W.210kgs