



LAMPE A FENTE

UT-5F Slit Lamp Microscope

Characteristics:

The optical lens are made of glass imported from Germany.
All the lens are multi-coated, mildewproof and anti-reflection.
The adaptor for Applanation Tonometer are universal.

Specifications:

Microscope Type	Converging Stereoscope
Magnification	three-step magnifications
Eye piece	12.5x
Total Magnification and Diameter of Visual Field(mm)	10x (Φ18) 、 16x (Φ14.5) 、 25x (Φ9)
Adjustable Range of Pupil Distance(mm)	55~75
Adjustable Diopter	-5D~+3D
Slit Width(mm)	0-9 adjustable
Slit Height(mm)	1~8 adjustable
Slit Angle	0~180°
Slit Inclination	5°、 10°、 15°、 20°、 4 ranges chosed
Light Spot Diameter (mm)	Φ9、 Φ8、 Φ5、 Φ3、 Φ2、 Φ1、 Φ0.2
Filter	heat absorption , grey ,red free , cobalt blue
Illumination Bulb	12V / 30W halogen bulb
Voltage/frequency	AC220V±10%/60Hz±1Hz. AC110V±10%/50Hz±1Hz
Power	58VA
Electrical Safety Standard	Conform to Standard IEC 601-1.Class1Type B
Whole Packing Volume	0.30m ³
Carton Number	2
Total Weight	51kg



Accessories:

Applanation Tonometer (R type and T type), Measuring Eyepiece