



.VIDEO PRINTER SONY UP-897MD

Thermal Head	Thin Film Thermal Head (1280 dots)
Gradations	True 256 gray balance
Resolution	325 DPI
Printed Picture Size	4"(W) x 12 5/8" (H)
Effective Pixels	EIA 1280 dot x 508 lines CCIR: 1280 dot x 608 lines
Memory	10 Frames
Input/Output	Video In: BNC, CCIR or EIA composite video signal, 1.0 Vp-p, 75 ohms/ high impedance selectable
Video Out:	BNC, Loop-through or E to E selectable
Print Speed	Approximately 2.0 seconds (High Speed Mode)
Remote Terminal	Stereo Mini-jack
Power Requirements	100V to 240V, 50/60 Hz, 1.5A to 0.8A (Worldwide)
Dimensions (WHD)	6 1/8 x 3 1/2 x 9 1/2 inches 154mm x 88mm x 240mm
Weight	Approximately 5 lbs. 12oz. 2.6kg
Safety Standards	UL60950-1 Listed + CSA C22.2 No. 60950-1 + UL60601-1 Listed + CSA No. 60601-1
Supplied Accessories	AC Power Cord Operation Manual (CD ROM) Thermal Head Cleaning Sheet BNC Video Connecting Cable UPP-110HG Media (1 Roll)
Media (Not Included)	UPP-110HG High Gloss Paper UPP-110HD High Density Paper UPP-110S Standard Paper

