



APPAREIL RX - TRANSPORTABLE SCF-10

The unit is shock-proof, single-focus and self-rectified X-ray unit.
It applies to medical fluoroscopy and radiography on mouth, carity bone , heart , skull as well as contraceptive ring.

Specifications:

Maximum rating:

Far fluoroscopy: 75 kVp 3-5mA
[uncontinuous use]

Far radiography: 75 kVp 10mA 6 sec,
(interval time: 3 minutes)

(Routine fluoroscopy condition: 65kVp 3mA)

Power supply: AC 220V,50HZ, 1kVA

Current: for fluoroscopy 2A(uncotinuuous) for radiography 6A(instantaneous)

Size of screen: 203mm x254mm

Range of movement :

up-down: 220mm,

left-right: 130mm

Tube turning angle:

vertical 360°,

horizontal 320°

Volume: 50cmx31cmx18cm

Weight: 25kg

